

SHRILL NOISES REDUCER

DENTAL

MAGi(C)
SONIC





SOUND PRESSURE DAMPER

MAGIC SONIC consists of 3 resonators which reduce the sound pressure in the frequency range of around 6000 Hz.

The sound environment and the perception of speech sounds remain unchanged for the comfort of dentists, hygienists, assistants and patients.

MAGIC SONIC does not create an occlusion effect (stuffy ear sensation).



PATENT

Made in France



The box contains :



**MAGIC
SONIC**

x 4
(2 pairs)

