

Finishing & Polishing

Sof-Lex™ Finishing Strips


Sof-Lex finishing and polishing strips are designed to produce a smooth restoration surface to help prevent interproximal plaque buildup for superior finishing and polishing results. Designed to complete the Sof-Lex system, Sof-Lex finishing strips are gapped at their centers for easy interproximal insertion and are precision-coated with an aluminum-oxide grit.

























Sof-Lex™ Finishing and Polishing Aluminium Oxide Strips

Product Information

Item Number

150 Centre Gapped Coarse/Medium	1954
100 Centre Gapped Coarse/Medium - Narrow	1954N
100 Centre Gapped Fine/Super Fine	1956



Abrasive Grades	Sof-Lex Moore's Mandrel Type	Sof-Lex	Sof-Lex Extra Thin (XT)	Sof-Lex Finishing Strip
Superfine	 1958SF	 1982SF  1981SF	 2382SF  2381SF	SF 
Fine	 1958F	 1982F  1981F	 2382F  2381F	F  1956
Medium	 1958M	 1982M  1981M	 2382M  2381M	M 
Coarse	 1958C	 1982C  1981C	 2382C  2381C	C  1954N 1954