

# EzSensor Classic™

## EASY DIAGNOSIS AND GREAT EXPERIENCE WITH EZSENSOR CLASSIC



EzSensor Classic consistently provides high resolution and clear images for accurate diagnosis and treatment planning.

### VARIOUS SIZES (1.0 / 1.5 / 2.0)

Select the right sized sensor from the three size options to increase patient comfort.

### ERGONOMIC DESIGN

The EzSensor Classic has a slim design with rounded corners for easy positioning to ensure patient comfort.



Square Corner



Rounded Corner

### DURABILITY

The EzSensor Classic's unique design makes itself extremely durable. The exterior is made of rugged aluminum and the interior is designed to absorb the external shock. Furthermore, a reinforced, flexible cable attachment protects the sensor from excessive bending.



### SPECIFICATIONS (EzSensor Classic : IOS-U10IB / IOS-U15IB / IOS-U20IB / IOS-U10VB / IOS-U15VB / IOS-U20VB)

Detector	CMOS
Pixel Size	29.6 μm
Theoretical Resolution	17 lp/mm
Dynamic Range	12 bit
Active Area (WxL)	Size 1.0 : 20X30 mm Size 1.5 : 24X33 mm Size 2.0 : 26X36 mm
Dimensions (WxLxT)	Size 1.0 : 25.4X36.8 mm (1.00"X1.45") Size 1.5 : 29.2X39.5 mm (1.14"X1.55") Size 2.0 : 31.3X42.9 mm (1.23"X1.69")
Thickness	4.8 mm (0.19")
Cable Length	2.7 m

\* The specifications are subject to change without prior notice.

### [ Intended use ]

Intra-Oral Sensor is intended to collect dental x-ray photons and convert them into electronic impulses that may be stored, viewed and manipulated for diagnostic use by dentists.